
Heraeus
KULZER

Technomat

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Pressure equipment for cold-curing resins used in materials testing

When fast-curing resins are allowed to set in a pressure unit, one obtains bubble- and pore-free samples. The physical and chemical properties of the resin however remain completely unchanged.

The Technomat is a pressure unit that was designed to suit the special technology of fast-curing resins. The pressure level of 2.0 bar was chosen to make allowance for the polymerization behaviour of the cold polymer Technovit, which is used for embedding metallographic specimens and for taking surface impressions in materials testing. In processing Technovit, one can achieve

- * clear, bubble-free, transparent specimens in the case of Technovit 4004,
- * pore-free embedded specimens in the case of Technovit 4000 and 4071,
- * bubble-free impressions in the case of Technovit 3040,
- * pore-free castings in the case of Technovit 2060

if curing is carried out in the Technomat.

The functional design of the unit makes it simple and problem-free to operate. When the valve lever is moved upwards, the pressure is built up to a level of 2.0 bar. The pressure gauge within the unit housing indicates the operating pressure. The pressure is released by moving the valve lever downwards. The Technomat is a space-saving, compact unit. The housing is made from macrolon resin and the pressure pot itself is of INOX steel.

Technical Data

Operating pressure	2,0 bar
Safety valve	2,8...3,3 bar
Pressure hose connection	3...10 bar
Dimensions	width 340 mm depth 340 mm height 255 mm
Weight	4,0 kg

Right of technical modification reserved.

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Hersteller/Manufacturer: Heraeus Kulzer GmbH
Heraeus Kulzer GmbH, Wehrheim